



Surgical Instruments

Crimson International is proudly the most sought after name in the field of Medical Surgical Instruments. Since its birth, the chief aim of our company has been to utilize all its energies, manpower and technology towards developing superior quality instruments, with special emphasis on them being technically precise as well as user friendly. In order to serve the interest and meet the requirement of our customers, we endeavor to best of our abilities to improve our existing line and also introduce new products.



19

Sterilization and Holloware



sizes in mm ø x h	capacity l	satin finish
165 x 60	0.75	DC19-075-16
175 x 70	1.00	DC19-075-17
220 x 80	2.00	DC19-075-22
225 x 90	2.50	DC19-075-23
240 x 110	3.00	DC19-075-24
255 x 115	4.00	DC19-075-25
275 x 115	4.50	DC19-075-27
290 x 125	6.00	DC19-075-29*
310 x 120	6.00	DC19-075-31
345 x 125	8.50	DC19-075-34
380 x 130	11.00	DC19-075-38
415 x 140	14.00	DC19-075-41



sizes in mm ø x h	capacity l	highly polished	satin finish
40 x 19	0.05		DC19-079-04
61 x 30	0.07		DC19-079-06
116 x 50	0.35	DC19-077-11	DC19-079-11
128 x 55	0.45	DC19-077-12	DC19-079-12
147 x 65	0.7	DC19-077-14	DC19-079-14
167 x 75	1.1	DC19-077-16	DC19-079-16
187 x 85	1.6	DC19-077-18	DC19-079-18
ø x h			
116 x 35	0.3	DC19-081-11	DC19-083-11
128 x 38	0.4	DC19-081-12	DC19-083-12
147 x 40	0.5	DC19-081-14	DC19-083-14
167 x 45	0.75	DC19-081-16	DC19-083-16
187 x 45	1.0	DC19-081-18	DC19-083-18



DC19-085-02 = 0.25l
 DC19-085-05 = 0.5l
 DC19-085-10 = 1.0l
 DC19-085-20 = 2.0l

Lid only
 DC19-091-10
 100 x 16 mm



Dish + Lid
 DC19-087-10
 100 x 20 mm

Dish only
 DC19-089-09
 96 x 20 mm



Sizes in mm
Ø x h

Jar without lid

Lid only

Lid only with triangle hole

30 x 90
50 x 100
50 x 130
50 x 175
75 x 100

DC19-093-30
DC19-093-50
DC19-093-51
DC19-093-52
DC19-093-75

-
DC19-095-50
DC19-095-50
DC19-095-50
DC19-095-75

-
DC19-097-50
DC19-097-50
DC19-097-50
DC19-097-75

75 x 130
75 x 175
100 x 175
100 x 200

DC19-093-76
DC19-093-77
DC19-093-98
DC19-093-99

DC19-095-75
DC19-095-75
DC19-095-99
DC19-095-99


DC19-097-75
DC19-097-75
DC19-097-99
DC19-097-99




DC19-099-23
235 x 190 x 40 mm




DC19-101-22
Ø 220 x 70 mm




sizes in mm l x b x h	lid with knob	lid without knob	lid with handle	lid with flush handle
165 x 085 x 35	DC19-103-16	DC19-105-16		
200 x 100 x 35	DC19-103-20	DC19-105-20		
200 x 100 x 50	DC19-103-21	DC19-105-21		
230 x 130 x 50	DC19-103-23	DC19-105-23		
260 x 150 x 50	DC19-103-26	DC19-105-26	DC19-107-26	DC19-109-26
270 x 235 x 50	DC19-103-27	DC19-105-27	DC19-107-27	DC19-109-27
300 x 200 x 50	DC19-103-30	DC19-105-30	DC19-107-30	DC19-109-30
325 x 375 x 30	DC19-103-32	DC19-105-32	DC19-107-32	DC19-109-32
325 x 275 x 50	DC19-103-33	DC19-105-33	DC19-107-33	DC19-109-33
420 x 080 x 40	DC19-103-42	DC19-105-42	DC19-107-42	DC19-109-42
420 x 175 x 20	DC19-103-43	DC19-105-43	DC19-107-43	DC19-109-43




sizes in mm Ø x h	lid with knob	with lid without knob
075 x 60	DC19-111-06	DC19-123-06
075 x 75	DC19-111-07	DC19-123-07
100 x 60	DC19-111-06	DC19-125-06
100 x 100	DC19-111-10	DC19-125-10
125 x 80	DC19-111-08	DC19-127-08
125 x 125	DC19-113-12	DC19-127-12
135 x 20	DC19-113-02	DC19-129-02
135 x 35	DC19-113-03	DC19-129-03
135 x 50	DC19-113-05	DC19-129-05
150 x 100	DC19-113-10	DC19-131-10
150 x 150	DC19-115-15	DC19-131-15
200 x 100	DC19-115-10	DC19-133-10
200 x 150	DC19-115-15	DC19-133-15
200 x 200	DC19-117-20	DC19-133-20
240 x 160	DC19-119-16	DC19-135-16
240 x 230	DC19-119-23	DC19-135-23
290 x 160	DC19-121-16	DC19-137-16
290 x 240	DC19-121-24	DC19-137-24




Phillip
DC19-139-01
ø 75 x 115 mm



Phillip
DC19-141-01
ø 75 x 100 mm
with Glass Insert



Phillip
DC19-141-05
ø 65 x 70 mm
Glass Insert only




ø x h
135 x 40 mm



Martin Sterilette 1
DC19-143-13 lid with knob
DC19-143-14 lid without knob



Martin Sterilette 2
DC19-145-13 lid with knob
DC19-145-14 lid without knob
DC19-147-13 Set of special filters = 5 pairs




l x b x h
200 x 100 x 40 mm



Martin Sterilette 3
DC19-149-20 lid with knob
DC19-149-21 lid without knob



Martin Sterilette 4
DC19-151-20 lid with knob
DC19-151-21 lid without knob
DC19-153-20 Set of special filters = 5 pairs



sizes in mm


l	b	h
50	25	4.5
65	42	8
75	25	4.5



DC19-155-50
DC19-155-65
DC19-155-75



DC19-157-50
DC19-157-65
DC19-157-75



complete Insert only

DC19-159-65	DC19-161-65
ø = 65 mm	ø = 62 mm
h = 15 mm	h = 13 mm



DC19-159-65



DC19-161-65



DC19-163-13
l x b x h
135 x 100 x 25 mm



DC19-165-35
35 mm high



DC19-165-50
50 mm high



DC19-165-60
60 mm high




All bottoms with filters

DC19-165-35 - Deponette with filter bottom, height 35 mm
DC19-165-50 - ditto, height 50 mm
DC19-165-60 - ditto, height 60 mm




Type 7
Disposal Container
Lid unperforated
Bottom unperforated

DC19-167-15



Type 8
Sterilizing box
Lid perforated
Bottom perforated
hinged Lid
Long-life textile-filter
Filter change by releasing
the Knurled nuts

DC19-169-07




DC19-171-12
125 x 46 x 15 mm



Schimmelbusch
DC19-173-08 - DC19-173-39

ø mm	h mm	
080	080	DC19-173-08
120	120	DC19-173-12
140	140	DC19-173-14
160	120	DC19-173-16
160	150	DC19-173-17
190	150	DC19-173-19
190	190	DC19-173-20
240	160	DC19-173-23
240	240	DC19-173-24
290	160	DC19-173-28
290	240	DC19-173-29
340	160	DC19-173-33
340	240	DC19-173-34
390	190	DC19-173-38
390	290	DC19-173-39




Sauerbruch
DC19-175-20
210 x 125 x 17 mm




½

Bunt
DC19-177-12
12 cm/4 ¾"



½

Mayo
DC19-179-14
14 cm/5 ½"



½

Mayo
DC19-181-14
14 cm/5 ½"